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#### Pressure transducers in the DPT-6000 range

Please read this carefully before use.

#### Application areas

The pressure transducers of the DPT-6000 range are used for invasive pressure measurement to measure the pressures of the human circulation system and convert these into electrical signals.

#### Product description

The pressure transducers in the DPT-6000 range are equipped with a highly sensitive and precise microchip for the transmission and reproduction of the physiological pressure curve. The converted pressure signals can be displayed and monitored on the patient monitor via a suitable interface cable.

In order to ensure unobstructed flow in the connected catheter or cannulae, some pressure transducer models are fitted with an integrated flush device. This guarantees a continuous and constant infusion of saline solution, whose flow rate is determined by the patient's blood pressure, the pressure of the pressure bag on the flush device and selected capillary. The continuous flow of the selected flush device (3ml/h, 30 ml/h) is triggered by a nominal pressure of 300mmHg.

#### Application

1 Prepare the measuring site by assembling the universal clamp and fixing the mounting plate. In the case of double and multiple monitoring sets, the mounting plates can be added on in building-block (modular) fashion.

2 Open the sterile packaging of the pressure transducer/monitoring set. Carefully insert the DPT-6000 pressure transducer into the guide of the mounting plate from the top until it stops.

3 For connection with the monitor, connect the splash-proof transducer plug into the interface cable. The interface cable is fitted with a special monitor plug depending on the type of monitor.

4 Turn on the monitor. Take note that electronic devices require at least 5 minutes to reach their operating temperature and work without „drifting“.

#### Filling the monitoring set

a) Close down the roller clamp on the infusion tubing. Then place the spike of the nonventilated drip chamber into the infusion bag filled with sterile saline solution using a light turning movement.

b) Make sure that filling takes place by force of gravity, i.e. that there is no pressurisation on the flush solution. If pressure is exerted before filling, this can lead to the formation of microbubbles which will prejudice the subsequent pressure measurement and thus lead to incorrect pressure values. If the flush flow rate is slow, after squeezing the flush device, then the infusion bag must be hung at a higher level.

c) In order to fill by the force of gravity, the infusion bag must be at a higher level than the pressure transducer.

d) Before the saline solution flows through the infusion set, the drip chamber must be at least 2/3 full (liquid level setting).

e) Close the patient pressure line with the stopcock on the pressure transducer (OFF arm of the tap points in the direction of the patient pressure line). Open the roller clamp and purge the pressure transducer including the side outlet of the zero-point stopcock by carefully squeezing the flush device wings. Remember that air rises. For this reason it should always be filled from a lower level.

f) Open the tap of the pressure transducer in the direction of the patient (OFF arm is pointing to the side outlet of the stopcock) and fill the patient line. Check the monitoring set for air bubbles. The pressure bag must only be pressurised with 300 mmHg of pressure when all of the air bubbles have been purged from the system. Continuous flow is now regulated by the flush device.

6 Remove the protective cap from the Luer lock connection of the patient line and allow a liquid drop to form on the vertically held Luer. Now connect firmly to the catheter or cannula connected to the patient.

#### Zero-point adjustment/Calibration

a) Position the pressure transducer in such a way that the side port of the integrated three-way stopcock is at the level of the organ to be monitored. If the pressure transducer is in a horizontal position the flow channel of the pressure transducer must be at the level of the organ to be examined.

b) Open the side port of the 3-way stopcock to atmosphere by turning the tap 90° upwards (OFF arm is pointing in the direction of the patient).

c) Carry out the zero-point calibration at the monitor in accordance with the monitor manufacturer's instructions. Then return the 3-way stopcock to its original position by turning it 90°.

8 The monitoring set is now ready for continuous pressure monitoring. If the system does not operate as described above, repeat the various steps. If the problem cannot be eliminated, replace the pressure monitoring system.

#### Technical data

Measurement range	-30 to +360 (mmHg)
Maximum permissible pressure	4000 mmHg
Sensitivity	5 µV/mmHg (±1%)
Temperature-related sensitivity	±0,1%/°C
Temperature-related zero-point drift	±0,1 mmHg/°C
Linearity and hysteresis error	±1%
Maximum zero-point error (offset)	≤25 mmHg
Ausgangsimpedanz	356 Ω ±10%
Betriebstemperatur-Bereich	+15 °C bis +40 °C
Versorgungsspannung	2 - 15 V oder ≈ (bis 5 kHz)
Patientenableiterstrom	≤10 µA
Light sensitivity at 100.000 lux	not verifiable
Monitor compatibility	compatible with all standard monitors BF/CF
Defibrillation resistance	5 discharges of 360 Joules at 50 Ω load in 5 minutes
Isolation of the electric components from the fluid channels	≥5000 Volt

#### Warnings

A The plug housing and the plug pins of the transducer cable must not come into contact with liquid or disinfecting agent, as this can cause corrosion of the pins.

B Before filling the system always check the electrical connection between the transducer cable and the interface cable.

C Take special care during the filling and flushing process that the electrical connection does not come into contact with fluids.

D If it comes into contact with liquid there is a danger that the pressure signals will be prejudiced which can lead to an incorrect reading of the pressure values or even a complete breakdown of the pressure transducer.

E Do not disconnect the pressure transducer from the interface cable socket by pulling on its cable. To disconnect the pressure transducer, pull the pin on the recessed grip out of the interface cable.

F The DPT-6000 series pressure transducers are designed for single use on one patient.

G Contents sterile unless package is opened or damaged.

H The set and the components must be disposed of with great care and in accordance with the relevant guidelines for contaminated medical waste.

I As is known from the literature, invasive pressure measurement can have subsequent complications: Over- and underdraining of the pressure signal, faulty zero-point calibration, thrombosis, accidental injection of medication, danger of haemorrhage when disconnecting, aerobolism, haematomas, infections and ruptures of the artery.

J If saline solution leaks from the side outlet of the transducer stopcock during the zero-point calibration (compare with points